

Product(s): Reveal® G4 Rapid HIV-1 Antibody Test
Subject: IQCP Risk Assessment and Quality Control Plan

Dear Customer,

MedMira has compiled this information to aid in the preparation of your laboratory’s IQCP Risk Assessment and Quality Control Plan for Reveal G4 Rapid HIV-1 Antibody Test (Reveal G4). This information is not intended to replace the Package Insert provided with the test. Information has been abstracted from the Reveal G4 Package Insert (IRIPYZIS0001EN, Rev 1/1) and the Reveal G4 HIV-1 Antibody Test Controls Package Insert (IRICPIS0001EN, Rev 0/1).

| Failure | Cause | Failure Type | Reveal G4 Risk Mitigation Features | Reveal G4 Package Insert Section | Potential Effect(s) of Failure | Laboratory Risk Mitigation | Severity | Frequency | Laboratory Documentation |
|----------------------|---|--------------|------------------------------------|---|--------------------------------|----------------------------|----------|-----------|--------------------------|
| Inappropriate sample | Substances interfering with sample flow into the Test Cartridge | Specimen | None | Limitations of the Test | | | | | |
| | Substances interfering with product formulation/biologics | Specimen | None | Performance Characteristics > Interfering Substances and Unrelated Medical Conditions | | | | | |
| | Non-validated sample type collected | Operator | None | Specimen Handling/ Collection and Use | | | | | |
| | Non-validated anticoagulant used during sample collection | Operator | None | Specimen Handling/ Collection and Use | | | | | |

| Failure | Cause | Failure Type | Reveal G4 Risk Mitigation Features | Reveal G4 Package Insert Section | Potential Effect(s) of Failure | Laboratory Risk Mitigation | Severity | Frequency | Laboratory Documentation |
|------------------------------|--|--------------|------------------------------------|--------------------------------------|--------------------------------|----------------------------|----------|-----------|--------------------------|
| Incorrect sampling procedure | Sample contaminated with other body fluids/materials during collection | Specimen | None | Precautions | | | | | |
| | Incorrect or non-validated sample type tested | User | None | Specimen Handling/Collection and Use | | | | | |
| | Unapproved patient category collected | User | None | Specimen Handling/Collection and Use | | | | | |
| | Inadequate sample volume collected | User | None | Specimen Handling/Collection and Use | | | | | |
| | Incorrect technique for sample collection | User | None | Specimen Handling/Collection and Use | | | | | |
| Incorrect sample handling | Sample contaminated after collection | User | None | Specimen Storage section of PI | | | | | |
| | Sample stored at incorrect temperature | Environment | None | Specimen Handling/Collection and Use | | | | | |
| | Sample stored for too long prior to testing | Environment | None | Specimen Handling/Collection and Use | | | | | |
| | Whole blood specimen frozen | User | None | Specimen Handling/Collection and Use | | | | | |
| | Frozen serum/plasma specimens mishandled during thawing | User | None | Specimen Handling/Collection and Use | | | | | |
| | Serum/plasma samples frozen and thawed more than twice | User | None | Specimen Handling/Collection and Use | | | | | |
| | Samples tested before equilibrating to room temperature | User | None | Specimen Handling/Collection and Use | | | | | |

| Failure | Cause | Failure Type | Reveal G4 Risk Mitigation Features | Reveal G4 Package Insert Section | Potential Effect(s) of Failure | Laboratory Risk Mitigation | Severity | Frequency | Laboratory Documentation |
|------------------|---|--------------|--|---|--------------------------------|----------------------------|----------|-----------|--------------------------|
| Operator failure | Package insert instructions not followed | User | None | Package Insert | | | | | |
| | Incorrect sample volume added to the Test Cartridge (Plasma) or Universal Buffer vial 1 (fingerstick Whole Blood) or Sample Tube (venipuncture Whole Blood) | User | None | Specimen Handling/ Collection and Use, Test Procedure | | | | | |
| | Patient sample used as external control material | User | None | Test Procedure | | | | | |
| | External control material substituted as patient sample | User | None | Test Procedure | | | | | |
| Reagent handling | Product damaged during shipment | Environment | None | Precautions > Handling Precautions | | | | | |
| | Storage temperature outside specified range (2-30°C) | Environment | None | Precautions > Handling Precautions | | | | | |
| | Used Mylar pouch from which Silica packet was missing | User | None | Precautions > Handling Precautions | | | | | |
| | Used Mylar pouch which was punctured or opened in advance of use | User | None | Precautions > Handling Precautions | | | | | |
| | Used Mylar pouch which has been used was damaged | User | None | Precautions > Handling Precautions | | | | | |
| | Reused Test Components (excluding Universal Buffer) | User | None | Precautions > Handling Precautions | | | | | |
| | Used Test Cartridge with contaminated or damaged immunoreactive membrane in Test Cartridge | User | None | Precautions > Handling Precautions | | | | | |
| | Used Expired reagents | User | Expiration printed on reagents and packaging | Precautions > Handling Precautions | | | | | |

| Failure | Cause | Failure Type | Reveal G4 Risk Mitigation Features | Reveal G4 Package Insert Section | Potential Effect(s) of Failure | Laboratory Risk Mitigation | Severity | Frequency | Laboratory Documentation |
|------------------------|--|--------------|---|------------------------------------|--------------------------------|----------------------------|----------|-----------|--------------------------|
| | Used reagents from different lots | User | Lot number printed on reagents and packaging | Precautions > Handling Precautions | | | | | |
| Reagent Test Procedure | General Reagent Failure | Reagent | Built-in Procedural and Reagent Control is designed to detect general reagent failures. | Quality Control | | | | | |
| | Sample not added to Test Cartridge | User | Built-in Procedural and Reagent Control is designed to detect procedural errors. | Test Procedure, Quality Control | | | | | |
| | Reagent deterioration | Reagent | None | Storage Instructions | | | | | |
| | Test Procedure not followed correctly | User | None | Test Procedure | | | | | |
| | Product used before equilibrating to room temperature | User | None | Test Procedure | | | | | |
| | Reagent contamination | User | Universal Buffer formulation includes preservative agent | Storage Instructions | | | | | |
| | Inappropriate/inadequate volume of Universal Buffer added to InstantGold cap | User | None | Test Procedure | | | | | |

| Failure | Cause | Failure Type | Reveal G4 Risk Mitigation Features | Reveal G4 Package Insert Section | Potential Effect(s) of Failure | Laboratory Risk Mitigation | Severity | Frequency | Laboratory Documentation |
|----------------------------|--|--------------|--|---|--------------------------------|----------------------------|----------|-----------|--------------------------|
| Test result Interpretation | Inadequate lighting used to read test results | User | None | Handling Precautions | | | | | |
| | Test Results not read immediately after completion of Test Procedure | User | None | Test Procedure, Test Result Interpretation | | | | | |
| | Misinterpretation of test results | User | None | Test Results and Interpretation of Results | | | | | |
| External Control failure | Shipping temperature outside control's specified range (2-30°C) | Environment | None | Handling and Storage Instructions (Test Control Package Insert) | | | | | |
| | Storage temperature outside specified range (2-30°C) | Environment | None | Handling and Storage Instructions (Test Control Package Insert) | | | | | |
| | External control materials tested before equilibrating to room temperature | User | None | Directions for Use > General Test Preparation (Test Control Package Insert) | | | | | |
| | Reactive and Non-Reactive Test Control misused | User | None | Directions for Use (Test Control Package Insert) | | | | | |
| | External control materials expired | User | Expiration printed on reagents and packaging | Directions for Use > General Test Preparation (Test Control Package Insert) | | | | | |

Should you have any questions regarding this Customer Advisory Notice or Reveal G4 please contact our Customer Support Team at G4@medmira.com or +1 877 633 6472.