

Eliminate the mess of cerumen removal during irrigation

# OtoClear® Ear Irrigation Tip

Safe, soft, comfortable 3-streamed-tip that prevents the risk of over-insertion



bionix®



MADE IN THE USA



## FEATURES

- Made with soft, flexible material for patient comfort
- Uses three divergent streams of water to safely loosen cerumen
- Directs water toward the ear canal wall instead of at the tympanic membrane
- Patented design prevents backsplash

## BENEFITS

- Compatible with a variety of devices
- Improves efficiency and patient safety
- Reduces patient discomfort
- Reduces messy clean up of cerumen removal during irrigation

## Bionix OtoClear Ear

**Irrigation Tip** is a soft, flared, 3-streamed-tip that improves patient safety and eliminates the mess of cerumen removal during irrigation.

It's compatible with a variety of delivery devices, meeting personal preferences of healthcare professionals.

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

 OTOCLEAR TIPS ARE SINGLE USE ONLY

 **Rx ONLY**



OTOCLEAR TIPS ARE MADE IN THE USA

Contact your Bionix Sales Rep, Distributor of Choice or call 1.800.551.7096 (US)  
International: Contact your Distributor or call +1 419.727.8421

©2022 Bionix LLC  
MKT-0041 | REV. E

INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

1670 Indian Wood Circle | Maumee, OH 43537  
United States | P: 800.551.7096 (+1 419.727.8421)  
F: 800.455.5678 (+1 419.727.8426) | www.bionix.com  
Email: bionix@bionix.com (international@bionix.com)

# OtoClear® Ear Irrigation Tip

## Product Offering



### OtoClear Ear Irrigation System\*

Recommended for high volume facilities

- A completely portable irrigation delivery device that can be used for procedures anywhere in or out of the office.
- Provides effective, pulsed irrigation.
- Higher volume than manual devices.
- Electronic pressure-control system with three settings.
- Proven safety and effectiveness.  
*(Children's Hospital of Pittsburgh Study)*
- Rechargeable battery.
- 278 ml self-contained water reservoir.

**#7280** - 1 Portable Electronic Irrigator, 40 OtoClear Tips, 1 Ear Basin

**#7245** - 1 Portable Electronic Irrigator, 5 OtoClear Tips, 1 Ear Basin

**\*PRODUCT AVAILABILITY VARIES BY COUNTRY**



### OtoClear SprayWash Kit

- Manual, intermittent lavage.
- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- 500ml bottle volume.

**#7290** - 1 SprayWash Bottle, 20 OtoClear Tips, 1 Ear Basin

**#7295** - 1 SprayWash Bottle, 5 OtoClear Tips, 1 Ear Basin



### OtoClear Ear Irrigation Tips

- Directs water to the ear canal wall, with three divergent streams instead of directly at the tympanic membrane.
- Flared design prevents backsplash and over-insertion.

**#7200** - OtoClear Ear Irrigation Tips (40 Tips per Box)

**#7202** - OtoClear Ear Irrigation Tips (20 Tips per Box)



### OtoClear Aquabot® Kit

- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- Provides continuous stream without manually pumping.
- 650ml bottle volume.

**#7270** - 1 Aquabot Bottle, 20 OtoClear Tips, 1 Ear Basin

**#7275** - 1 Aquabot Bottle, 5 OtoClear Tips, 1 Ear Basin

Aquabot is registered under Lunatech®, used under permission of Lunatech, Inc. All rights reserved.

Contact your Bionix Sales Rep, Distributor of Choice or call 1.800.551.7096 (US)  
International: Contact your Distributor or call +1 419.727.8421