

Lighted Ear Curette™ with Magnification

Bionix® Lighted Ear Curettes combine illumination and magnification to help healthcare providers take less time to see and safely remove cerumen from the ear canal.



PROBLEM

Ear curettage requires three hands. One hand to brace the patients head and straighten the ear canal, one to maneuver the curette and one to position a light to avoid blindly digging in the ear.

SOLUTION

The Lighted Ear Curette with Magnification brilliantly illuminates the ear canal eliminating the practice of blindly curetting the sensitive area.

RESULTS

Complete visualization shortens procedure time and increases efficacy.

“Bionix Lighted Ear Curettes have really changed our practice. We can quickly and painlessly clean out ears. We love them!” - Nancy Van Keuls, M.D.

POTENTIAL REIMBURSEMENT (U.S. ONLY)

CPT Code/Description

69210 - For potential reimbursement of impacted cerumen.



Obelis s.a • Boulevard Général Wahis 53 • 1030 Brussels, BELGIUM
Tel: +(32) 2.732.59.54 • Fax: +(32) 2.732.60.03 • E-Mail: mail@obelis.net
Bionix LLC • 1670 Indian Wood Circle, Maumee, Ohio 43537, United States
Phone: +1 419.727.8421 • Fax: +1 800.455.5678 • www.Bionix.com



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • www.bionix.com

©2021 Bionix LLC | MKT-0034 | Rev. D

Lighted Ear Curette™ with Magnification

PRODUCT OFFERING



Lighted AngleLoop® | #2210

Angled 4mm diameter loop, similar to Red AngleLoop® Safe Ear Curette™.



Lighted MicroLoop® | #2220

3mm diameter loop, similar to Green MicroLoop® Safe Ear Curette™.



Lighted WaveCurette™ | #2230

New wave-like 2.5mm wide tip.



Lighted FlexLoop® | #2240

4mm wide tip, similar to White FlexLoop® Safe Ear Curette™.



Lighted CeraSpoon® | #2250

4mm wide tip, similar to Yellow CeraSpoon® Safe Ear Curette™.



Lighted InfantScoop® | #2260

2mm wide tip, similar to Blue InfantScoop® Safe Ear Curette™.



Lighted VersaLoop® | #2270

3mm wide tip, similar to Purple VersaLoop® Safe Ear Curette™.



Lighted CeraPik® | #2280

2.5mm concave tip bent to 45°. Great for hard, impacted cerumen and foreign bodies.

STANDARD PACKS

(50 per box + 1 Light Source + 1 Magnification lens)

- #2210 - Lighted AngleLoop®
- #2220 - Lighted MicroLoop®
- #2230 - Lighted WaveCurette™
- #2240 - Lighted FlexLoop®
- #2250 - Lighted CeraSpoon®
- #2260 - Lighted InfantScoop®
- #2270 - Lighted VersaLoop®
- #2280 - Lighted CeraPik®
- #2245 - Variety Pack (Includes all tip styles except CeraPik®)

CLINIC PACKS

(200 per box + 4 Light Sources + 4 Magnification lenses)

- #2201 - Lighted AngleLoop® Clinic Pack
- #2202 - Lighted MicroLoop® Clinic Pack
- #2203 - Lighted WaveCurette™ Clinic Pack
- #2204 - Lighted FlexLoop® Clinic Pack
- #2205 - Lighted CeraSpoon® Clinic Pack
- #2206 - Lighted InfantScoop® Clinic Pack
- #2207 - Lighted VersaLoop® Clinic Pack
- #2208 - Lighted CeraPik®
- #2209 - Variety Pack (Includes all tip styles except CeraPik®)

“The variety of sizes and tip shapes makes ear wax removal easier than ever before.” - Carl A. Moritz, Jr., MD

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • www.bionix.com