

DOA Test Menu for USA

Unit: ng/mL

Drug	Abbreviation	Urine - Regular		Urine - Special			Oral Fluid		
		CLIA Waived	Forensic Use	CLIA Waived	510 (k) Cleared	Forensic Use	E&I (Employment & Insurance)	Forensic Regular	Forensic Special
Amphetamine	AMP	1000		500	300		50	50	40
Barbiturates	BAR	300					20	60	
Benzodiazepines	BZO	300			200	100 ♦		30	
Buprenorphine	BUP	10				5		5	
Cocaine	COC	300		150	100		20	20	
Cotinine	COT		200		200 ★★★				
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	EDDP	300				300			
Ethyl Glucuronide	EtG		500						
Fentanyl	FTY ★★★		20					30	
Gabapentin	GAB		2000 ^{NEW}						
Hydromorphone	HMO		300 ^{NEW}						
Ketamine	KET					1000			
Kratom	KRA		300 ^{NEW}			100 ^{NEW}			
Lysergic Acid Diethylamide	LSD		20 ^{NEW}						
Marijuana	THC	50 ★			40 ★	15 ^{NEW} & 25 ★	40★★	25 ★	
Methadone	MTD	300			200		30	30	
Methamphetamine	MET/mAMP	1000		500	300 ♦		50	50	
Methaqualone	MQL		300						
Methylenedioxymethamphetamine	MDMA	500					50	100	
Monoacetylmorphine	6-MAM		10						
Morphine	MOP/OPI300	300			100				
Opiate	OPI	2000					40	40	15
Oxycodone	OXY	100					20	20	
Phencyclidine	PCP	25					10	10	
Pregabalin	PGB		500 ^{NEW}						
Propoxyphene	PPX	300							
Synthetic Cannabinoid	K2		50						
Tramadol	TRA		200			100 & 1000 ♦			
Tricyclic Antidepressants	TCA	1000							
Alcohol	ACL (Saliva)							0.02 g/dL	
Alcohol	EtOH (Urine)					0.04 g/dL			

Adulteration Parameters: Creatinine (CR), Glutaraldehyde (GL), Specific Gravity (S.G.), pH, Oxidant (OX) & Nitrite (NI)

★ 11-nor-Δ⁹-THC-9-COOH

★★ Δ⁹-THC

★★★ Class I Exempt; Prescription Use Only

★★★★ Norfentanyl

◆ Please consult with Wondfo team before ordering